

URINE TESTING ALIQUOT REQUIREMENTS


Aliquot Container

URTS



Aliquot Transfer Tube
Lawson order number 004411

URCP



Mother Cup
Lawson order number 035712

URCN



Urine Conical
Lawson order number
Conical 063240 Caps 061974

UR60



Referrals
Order through Client supply

Test Codes or Workstations Associated with the Container

BKVUR MABSP COTN TP24 CRE24 CRCLR
TPURI/TPUR/PEPU* TPCUR CRUR TXPRO/TXSCR
WORKSTATION CHC01

URCUL
NONGY

PEPU/PEPFX/URFX* URSCR UAPR URSED
TXPRO/TXSCR* HCGU Workstation IMC05

WORKSTATION RLC26
WORKSTATION RLC49

Aliquot Container

VY



Viral Yellow

UR60-LP



Referrals – light protect

URTM



RLC30 – Med Tox

URG



RLC70– Green Top

Test Codes or Workstations Associated with the Container

CT/GC VAGPN
Order through Client Supply

WORKSTATION RLC26/60 when light
protection is called for
Order through client supply

WORKSTATION RLC30
Order through Client supply

WORKSTATION RLC70
(All toxicology quant levels)
Order through Client Supply

***TXPRO/TXSCR need 1 each of the white aliquot (URTS) and the Yellow Conical (URCN) (remember to only fill the URCN ¾ full)**

***PEPU/PEPFX WORKSTATION IMC05 needs to be aliquoted into BOTH the white Aliquot (URTS) and the yellow Conical (URCN)**